

# Isofield Alpha Kits

Sterile coverall with mandarin collar, hood & overboots



### PRODUCT INFO

### MANDARIN COLLAR COVERALL

GMP Grade A & B and ISO Class 4 & 5 cleanroom compatible Two-layer mandarin collar provides secure fit at neck area Zip flap with double-sided tape prevents ingress and egress of particles

Low-linting, anti-static fabric suitable for chemo handling Elasticated thumb loops

Back elastic to prevent bellows effect

### HOOD

Latex-free enclosed elastic face opening for comfortable fit Extra-long yoke for added protection Three panel hood design to minimize skin exposure

### **OVERBOOTS**

Extra-long knee height Two-ties for added security Anti-slip PU sole



All three items are aseptically folded

### **OUALITY ASSURANCE**

Manufactured in a facility operating under ISO 9001:2015 quality management system

Processed in a NEBB certified ISO Class 5 Cleanroom

### **APPLICATIONS**

Pharmaceutical product preparation Cleanroom cleaning and preparation Aseptic filling and labelling Medical device manufacturing Blending operations for biotechnologies Handling of active materials Maintenance of equipment and instruments











Isofield Alpha Kits



# MATERIAL PERFORMANCE PARAMETERS

2					
	DESCRIPTION	UNITS	METHOD	SPECIFICATION	TYPICAL
	Basic Weight	g/m²		±2	60
	Machine Direction Strength	N/5cm	ISO 9073-1-1989	≥50	69
	Machine Direction Elongation	%		≥100	109
	Cross Direction Strength	N/5cm		≥25	43
	Cross Direction Elongation	%		≥100	111
	Hydrostatic Head	mmH <sub>2</sub> 0	ISO 811-1981	≥100	213
• 1	Air Permeability	I/m²/S	GB/T5453	≥	2000
	Surface Resistivity	Ω	ISO 9073-1-1989	/	10-7
	Overboots Slip Resistance	COF	ASTM D1894	≥0.5	0.6 (Static) 0.5 (Dynamic)

CLEANLINESS LEVELS	UNITS	METHOD	SPECIFICATION	TYPICAL
Particle Shedding	≥ 0.3µm (counts/ min)	IEST-RP-CC003.4	2000	1000

MEINUD	CHEMOTHERAPY DRUGS TESTED	MEAN BREAKTHROUGH TIME (MBT), MINUTES BREAKTHROUGH OF THE TEST CHEMICAL IS DEEMED TO HAVE OCCURRED WHEN THE PERMEATION RATE HAS REACHED 0.1 µg/cm² /min*
	Carmustine	30 mins
ES	Cisplatin	> 480 mins
45 I M F / 39-20	Cyclophosphamide	> 480 mins
	Doxorubicin HCI	> 480 mins
	Etoposide	> 480 mins
	Fluorouracil	> 480 mins
A	Methotrexate	> 480 mins
	Paclitaxel	> 480 mins
	Thiotepa	60 mins

<sup>\*</sup>Results achieved under controlled laboratory conditions, by accredited ISO 17025 third party testing laboratory

### **PACKAGING**

1 set (hood/coverall/boots) per sealed inner bag 1 inner bag per sealed outer bag 40 outer bags per carton liner with twist tie (S-XXL) 35 sets per carton (3XL)

### **SHELF LIFE**

Three (3) years from date of manufacture

### **STORAGE**

Store in a cool, dry place ( $<40^{\circ}$ C) away from direct sunlight

Do not expose opened cartons to prolonged direct fluorescent light

**COUNTRY OF ORIGIN: CHINA** 

HS CODE: 6210109000

### PRODUCT CODE

## **REF** 80551

 SMALL
 80551S

 MEDIUM
 80551M

 LARGE
 80551L

 XLARGE
 80551XL

 XXLARGE
 80551XXL

 3XLARGE
 80551XXXL

Note: Overboots & hoods are universal size

Isofield Alpha Kits Page 3 of 6

# Isofield Beta Kits

Sterile coverall with mandarin collar, hood with integrated facemask & overboots



### PRODUCT INFO

### MANDARIN COLLAR COVERALL

GMP Grade A & B, ISO Class 4 & 5 cleanroom compatible Two-layer mandarin collar provides secure fit at neck area Zip flap with double-sided tape prevents ingress and egress of particles

Low-linting, anti-static fabric suitable for chemo handling Elasticated thumb loops

Back elastic to prevent bellows effect

### HOOD WITH INTEGRATED FACEMASK

Integrated mask with ties for quick and easy donning Latex-free enclosed elastic face opening for comfortable fit Extra-long yoke for added protection

Three panel hood design to minimize skin exposure

### **OVERBOOTS**

Extra-long knee height Two-ties for added security Anti-slip PU sole



All three items are aseptically folded

### **OUALITY ASSURANCE**

Manufactured in a facility operating under ISO 9001:2015 quality management system

Processed in a NEBB certified ISO Class 5 Cleanroom

### **APPLICATIONS**

Pharmaceutical product preparation Cleanroom cleaning and preparation Aseptic filling and labelling Medical device manufacturing Blending operations for biotechnologies Handling of active materials Maintenance of equipment and instruments











Isofield Beta Kits



# MATERIAL PERFORMANCE PARAMETERS

	DESCRIPTION	UNITS	METHOD	SPECIFICATION	TYPICAL
	Basic Weight	g/m²		±2	60
	Machine Direction Strength	N/5cm	ISO 9073-1-1989	≥50	69
	Machine Direction Elongation	%		≥100	109
	Cross Direction Strength	N/5cm		≥25	43
	Cross Direction Elongation	%		≥100	111
	Hydrostatic Head	mmH <sub>2</sub> 0	ISO 811-1981	≥100	213
	Air Permeability	I/m²/S	GB/T5453	≥	2000
	Surface Resistivity	Ω	ISO 9073-1-1989	/	10-7
	Overboots Slip Resistance	COF	ASTM D1894	≥0.5	0.6 (Static) 0.5 (Dynamic)

CLEANLINESS LEVELS	UNITS	METHOD	SPECIFICATION	TYPICAL
Particle Shedding	≥ 0.3µm (counts/ min)	IEST-RP-CC003.4	2000	1000

MEIHOD	CHEMOTHERAPY DRUGS TESTED	MEAN BREAKTHROUGH TIME (MBT), MINUTES BREAKTHROUGH OF THE TEST CHEMICAL IS DEEMED TO HAVE OCCURRED WHEN THE PERMEATION RATE HAS REACHED 0.1 µg/cm² /min*
	Carmustine	30 mins
2	Cisplatin	> 480 mins
<b>07-</b>	Cyclophosphamide	> 480 mins
IM F/39	Doxorubicin HCI	> 480 mins
	Etoposide	> 480 mins
	Fluorouracil	> 480 mins
A	Methotrexate	> 480 mins
	Paclitaxel	> 480 mins
	Thiotepa	60 mins

<sup>\*</sup>Results achieved under controlled laboratory conditions, by accredited ISO 17025 third party testing laboratory

### **PACKAGING**

1 set (hood/coverall/boots) per sealed inner bag, 1 inner bag per sealed outer bag,

40 outer bags per carton liner with twist tie (S-XXL), 35 sets per carton (3XL)

### **SHELF LIFE**

Three (3) years from date of manufacture

### **STORAGE**

Store in a cool, dry place ( $<40^{\circ}$ C) away from direct sunlight

Do not expose opened cartons to prolonged direct fluorescent light

**COUNTRY OF ORIGIN: CHINA** 

HS CODE: 6210109000

### PRODUCT CODE

## **REF** 80451

 SMALL
 80451S

 MEDIUM
 80451M

 LARGE
 80451L

 XLARGE
 80451XL

 XXLARGE
 80451XXL

 3XLARGE
 80451XXXL

Note: Overboots and hood are universal size

Isofield Beta Kits Page 5 of 6

# **Isofield Cleanroom Kits**





Isofield Cleanroom Kits